

Concrete Cement Rapid Accelerated Water Curing Tank

- Model NO. : Type B
- Temperature Range : Ambient to 100°C
- Transport Package : Standard Imports Packages
- Trademark : C-TECH
- Origin : China Mainland
- Manufacturer Company: shaoxing roan measuring & control technology co. ltd
- Importer and Bangladesh supplier : Rahim Trading Corporation -RTCDHAKA

Product Description : STANDARD: ASTM C684 This special curing tank has been designed for hot water curing in accelerated strength concrete. The interior is made from stainless steel, it is fitted with electronic programmer capable of controlling different test cycles with a choice of thermal gradients and curing time to a defined temperature value for a complete automatic curing cycle.

Temperature range	Ambient to 100°C
Inside dimension(W×D×H)	800×580×400mm
Capacity of the tank	2 samples 150×150mm cube mould 3 samples 100×100mm cube mould
Temperature Control Type	Automatic
Max Control Time	9999 min
Heat Power	10 kW
Power	AC380V, 50Hz

