

# Density Specific Gravity Balance

## Parameters

|                 |              |
|-----------------|--------------|
| Model           | HB-2000J     |
| Capacity        | 2000g        |
| Accuracy        | 0.01g        |
| Linearity Error | ±0.3g        |
| Pan size        | 195*195mm    |
| Dimension       | 290*205*70mm |
| Container Size  | Φ200*210mm   |
| Bracket Size    | Φ135*150mm   |
| Working Temp    | 5-40         |
| Power           | 220V         |
| Carton Size     | 440*400*460  |
| N/G.W           | 6kg/8kg      |



## Features :

- LCD display
- Aluminum base
- ABS body surrounding
- Stainless steel weighing pan
- Overload alarm/Error indication
- Leveling bubble/ Adjustable leveling feet
- RS3232 interface (optional)
- Printer(optional)
- Suitable for schools, research institute

